

## Generation of gRNAs for multiplexed targeting

This protocol is for synthesizing multiple guide RNAs for G0 phenotypic screening. It is adapted from PMID: 27809318 and uses Phusion polymerase, the Zymo DNA Clean and Concentrator-5 kit, and the Zymo RNA Clean and Concentrator-5 kit.

Note: Please ensure your top strand oligos are compatible with the bottom strand ultramer. The bottom strand ultramer sequence is:

AAAAGCACCGACTCGGTGCCACTTTTTCAAGTTGATAACGGACTAGCCTTATTTTAACTTGCTATTTCTAGCTCTAAAAC

### Assemble oligos

1. Prepare working stocks of top oligos with water to a concentration of 10  $\mu\text{M}$ .
2. Mix equimolar concentration of top oligos (10  $\mu\text{M}$  each). Up to 12 can be used in one reaction. E.g. mix 10  $\mu\text{l}$  of each oligo together. Vortex and brief centrifuge.
3. Set up the following oligo assembly reaction in a PCR strip:

2 x Phusion buffer +dNTPs	25 $\mu\text{l}$	1x
Top strand oligo mix (10 $\mu\text{M}$ )	1 $\mu\text{l}$	0.2 $\mu\text{M}$
Bottom strand ultramer (10 $\mu\text{M}$ )	1 $\mu\text{l}$	0.2 $\mu\text{M}$
Phusion DNA polymerase	0.5 $\mu\text{l}$	
Nuclease-free water	22.5 $\mu\text{l}$	
<b>Final</b>	<b>50 <math>\mu\text{l}</math></b>	

Note: the Bottom strand ultramer is found in the communal gRNA reagents box.

If synthesizing several, first make up a mastermix containing everything except the top strand oligo mix.

4. Anneal and extend using the thermocycler:

98°C 2 min  
 50°C 10 min  
 72°C 10 min  
 4°C hold

Note: this program is in the **General** folder, named **GRNA\_ASS**

### Purify assembled oligos (DNA Clean & Concentrator-5 kit, Zymo Research)

5. Transfer assembled oligo from PCR strip to 1.5 ml tube and add 250  $\mu\text{l}$  of DNA binding buffer.
6. Mix briefly by vortexing then spin down briefly to collect contents in bottom of tube.
7. Transfer mixture to a Zymo-Spin Column in a collection tube.
8. Centrifuge in desktop centrifuge at max speed for 30 seconds. Discard flow through.

9. Add 200  $\mu$ l of DNA Wash Buffer to the column (check ethanol has been added to the wash buffer). Centrifuge at max speed for 30 seconds and discard flow through.
10. Add another 200  $\mu$ l of DNA Wash Buffer, centrifuge for 30 seconds and discard flow through.
11. Transfer column to a fresh 1.5 ml tube (RNase-free); check there is no residual ethanol by dabbing the bottom of the column on your glove. If there is, spin again.
12. Add 6  $\mu$ l of RNase-free water directly to the column matrix. Incubate at room temperature for one minute.
13. Centrifuge for 30 seconds to elute the DNA.

### Quality Control: agarose gel

14. Run 1  $\mu$ l of purified assembled oligos on a 2% (wt/vol) agarose gel alongside 100 bp or 1 kb plus ladder. Correctly assembled oligos appear as 117-120 bp bands. This is an optional step and not necessary for experienced users.

### *In vitro* RNA synthesis

*From this point forward, use RNase-free solutions, tips, and technique. Clean work space with RNaseZAP, use fresh gloves etc.*

15. Using NEB T7 High Yield kit set up an NTP Buffer Mastermix containing buffer and each of the NTPs:

10 x buffer	1 $\mu$ l	2x
NTPs (100 mM)	1 $\mu$ l	20 mM each
<b>Final</b>	<b>5 <math>\mu</math>l</b>	

16. Set up T7 reaction mix in a PCR strip:

2 x NTP buffer mix	5 $\mu$ l	1x
Assembled oligo	4 $\mu$ l	
T7 RNA polymerase	1 $\mu$ l	
<b>Final</b>	<b>10 <math>\mu</math>l</b>	

17. Incubate for 4 hours at 37°C.

17. Add 0.5-1  $\mu$ L DNase

18. Incubate at 37°C for 15 mins, then bring volume to 50  $\mu$ L with nuclease-free water

### Purify sgRNA (RNA Clean & Concentrator-5 Kit, Zymo Research)

18. Transfer the contents of the in vitro RNA synthesis to a 1.5 ml RNase-free microfuge tube.
19. Add 2 volumes of RNA binding buffer (e.g., 100  $\mu$ L into a 50  $\mu$ L rxn)

20. Add an equal volume of 200 proof EtOH (e.g., 150  $\mu$ L). Mix by agitating the solution. Briefly spin to collect contents.
21. Transfer to a Zymo-spin column and collection tube, centrifuge for 30 s, discard flow-through
22. Add 400  $\mu$ L of RNA prep buffer, spin 30 s, discard flow-through
23. Add 700  $\mu$ L RNA wash buffer (check ethanol has been added), spin 30 s, discard flow-through
24. Add 400  $\mu$ L RNA wash buffer, spin 30 s, discard flow through
25. Spin dry for 2 min
26. Transfer column to a fresh RNase-free microfuge tube (check there is no residual ethanol), and elute in 15  $\mu$ L DNase/RNase free H<sub>2</sub>O by 30 s centrifugation.
  
27. Separate into aliquots (e.g., 5 aliquots 3  $\mu$ L each), flash freeze all but 1 or 2 in liquid N<sub>2</sub> or dry ice and store at -80°C.

### Quality Control: agarose gel and nanodrop

28. Run 1  $\mu$ L (or 1 aliquot) of purified RNA on a 2% (wt/vol) agarose gel alongside 100 bp or 1 kb plus ladder. Correctly assembled RNA appears as two 117-120 bp bands.
  
29. Nanodrop 1  $\mu$ L of purified RNA to obtain concentration. 260/280 should be  $\geq 2$ .